

Veterinary Anesthesia

COMPACT & STRONG

Moduflex Elite

Veterinary Anesthesia Machine

The Dispomed Moduflex Elite Veterinary Anesthesia Machine is ideal for rebreathing and small animal applications. Its compact size and base shape, combined with its large castors, provide excellent mobility & stability.

FEATURES

- 0.1 to 4 LPM O₂ flowmeter (up to 2 flowmeters can be installed).
- Rotary float and large scale for easy readability.
- Recessed flush valve prevents accidental activation.
- Quick connect fresh gas outlet for easy and secure switching between breathing circuits.
- SELECTATEC™ compatibility manifold for easy and quick installation of vaporizer
- Condensation trap to prevent water build-up in sodalime. Easy and efficient gas sealing.
- Momentary closed feature with pop-off valve.
- Standard CO₂ absorber
- Convenient anesthesia monitor shelf.
- 5 year warranty.



Premier XP Portable Anesthesia Machine

Standard Included Accessories

- Y Patient Circuit
- One 2 liter Rebreathing Bag
- One 3 liter Rebreathing Bag
- Fresh gas tube
- krytox grease
- Evacuation tube (150cm/5')

Product dimensions

| | |
|--------|-------------------|
| Length | 53.34 cm (21") |
| Width | 49.53 cm (19.50") |
| Height | 134.62 cm (53") |
| Weight | 25.90 kg (57 lbs) |

Shipping box dimensions

| | |
|--------|-------------------|
| Length | 71.12 cm (28") |
| Width | 68.58 cm (27") |
| Height | 53.34 cm (21") |
| Weight | 32.72 kg (72 lbs) |

Shown with optional vaporizer



Optional Accessories



Dual Cylinder Holder



Sodalime Usage



Circuit Pressure Alarm



Bain Circuit Adapter



Scavenger Interface



Maintenance Kit