

FOLDING WALL EXAM TABLES

INFINIUM  VET

PREMIUM SERIES

Constructed of 16 ga. Type 304 stainless steel 1 1/2" Sq. stainless square tubular legs, with 12 ga. wall supports. Easily folds up and out of the way when not in use creating extra space in your exam room for other various task. Top measures 44x22 and has a slightly raised marine edge to help prevent spills or accidents from going onto the floor area.

Standard SERIES

This wall mounted exam table is constructed of heavy-duty 18 gauge 304 stainless steel. The leg braces are constructed of 1-1/2" square stainless steel tubing and have extra heavy-duty 12 gauge stainless steel mounting brackets. This exam table includes a 2" high splash to help protect walls. Easily folds up and out of the way when not in use, creating extra space for other tasks.



INFINIUM[®]

Infinium Veterinary
12151 62nd St North #5
Largo, FL 33773 • USA

 1(866) 846-7517
 infiniumvet.com
 sales@infiniumvet.com