

Exam Lift Tables

LowMax Exam ET Electric Lift Table - 44" Top

PetLift Electric Veterinary Exam Table | Extra Low Height is a great exam table with an electric operation for graceful raising and lowering of the patient. This high-quality table boasts a smooth electric operation that lowers to 13.5" and raises to 41" high.

Large dogs walk right on. For added versatility, we have added a convenient duplex electric outlet on the base, for easy plug-in of electrically operated implements. This table has a sleek, slimline look and a powerful lift capacity. Perfect for small exam rooms.

FEATURES

- Product Weight : 260.0 LB
- Exam tabletop fabricated with 16 gauge, high quality 304 grade stainless steel.
- Top measures 22"W x 58"L with raised edge to prevent dripping.
- The base is double powder coated and will stand up to cleaning and disinfecting solutions.
- Foot operated pedals on both sides, raises from 13.5" to 41" high.
- Smooth lowering and lifting function.
- Capacity up to 300 lbs.
- A stationary table with plug-in electric operation.
- Duplex electric outlet on base.
- Adjustable leveling feet on base.
- Table weight is 295 lb.
- Made in the USA.



PN: 000.946.0015

POWER LIFT WITH FOOT CONTROL



EASY ACCESS: IMPRESSIVE LOW HEIGHT ADJUSTMENT



PN: 000.946.0080

OPTIONAL: NONSTICK TABLE TOPPER



Stainless LowMax Lift Table Specifications



VetLine

LowMax Electric Exam & Lift Table

Base Info

Lowers to: 13 1/2"

Raises to: 41"

Capacity: 300 lbs.

Power: Electric

Material: 7 ga.

Finish: Powder Coated

Product Dims & Weight

Width: 22"

Length: 44"

Weight: 200 lbs.

Top Specs

Width: 22"

Length: 44"

Material: 16 ga ss

Type: Fixed

Shape: Rectangle

Shipping Dims & Weights

Width: 40"

Length: 64"

Height: 30"

Weight: 295 lbs.

Vetline Series Comparison Chart

| PETLIFT LINE: | VetLine Series | VetLine Series | VetLine Series | VetLine Series |
|------------------------------------|----------------------|----------------------|----------------------|----------------|
| PRODUCT: | The Mobile Exam | Economy Exam HE | LowMax Exam ET | Classic Exam |
| BASE INFO | | | | |
| Lowers to | 17 | 19 | 13 1/2 | 27 |
| Raises to | 44 | 44 | 41 | 37 |
| Weight Capacity | 300 | 250 | 300 | 200 |
| Power | Battery | Hydraulic | Electric | Hydraulic |
| BASE SPECS | | | | |
| Base Material | 7 ga. | 7 ga. | 7 ga. | 7 ga. |
| Base Finish | Powder Coated | Powder Coated | Powder Coated | Laminate |
| Standard Color | Black or Silver Vein | Black or Silver Vein | Black or Silver Vein | Natural Almond |
| TOP SPECS | | | | |
| Left Hand(LH)/Right Hand(RH) | N/A | N/A | N/A | N/A |
| Width | 22 | 22 | 22 | 22 |
| Length | 58 | 44 | 44/58 | 44 |
| Shape | Rectangle | Rectangle | Rectangle | Rectangle |
| Top Material | 16 ga ss | 16 ga ss | 16 ga ss | 16 ga ss |
| Top Type | Fixed | Fixed | Fixed | Fixed |
| OVERALL DIMENSIONS | | | | |
| Width | 22 | 22 | 22 | 22 |
| Length | 58 | 44 | 44/58 | 40 |
| CABINET SPECS | | | | |
| No. of Drawers | N/A | N/A | N/A | 1 |
| No. of Doors | N/A | N/A | N/A | N/A |
| No. of Adjust Shelves | N/A | N/A | N/A | N/A |
| Miscellaneous | N/A | N/A | N/A | Rotating Base |
| ACCESSORIES INCLUDED | | | | |
| Duplex Electrical Outlet on Base | No | No | Yes | No |
| Corner Bump Guards A19-BG18 | Yes | No | No | No |
| Tie Downs A18-TD53 | Yes | No | No | No |
| Adjust Floor Levelers | No | Yes | Yes | No |
| Locking Casters A20-100-45 | Yes | No | No | No |
| Digital Scale 400 lbs. A16-PS500 | No | No | No | No |
| ACCESSORIES AVAILABLE | | | | |
| Duplex Electrical Outlet on Base | No | No | Incl. | No |
| Corner Bump Guards A19-BG18 | Incl. | No | No | Yes |
| Tie Downs A18-TD53 | Incl. | Yes | Yes | Yes |
| Adjust Floor Levelers A22-FL54 | No | Incl. | Incl. | No |
| Locking Casters A20-100-45 | Incl. | Yes | No | No |
| Waste Push Thru Flap A-20-100-45 | No | No | No | No |
| Table Mats | Yes | Yes | Yes | Yes |
| COLOR OPTIONS | | | | |
| Silver Vein | TEZ-9010BT-SV | TEZ-9025H-SV | TEZ-9095E-SV | N/A |
| Black | TEZ-9010BT-BLK | TEZ-9025H-BLK | TEZ-9095E-BLK | TBP-1010H-BLK |
| Natural Almond | N/A | N/A | N/A | TBP-1010H-NA |
| Designer White | N/A | N/A | N/A | TBP-1010H-DW |
| Fashion Gray | N/A | N/A | N/A | TBP-1010H-FG |
| Wood Grain: Kensington Maple | N/A | N/A | N/A | TBP-1010H-WGKM |
| Wood Grain: Wild Cherry | N/A | N/A | N/A | TBP-1010H-WGWC |
| SHIPPING DIMS & WEIGHTS | | | | |
| Width | 40 | 40 | 40 | 40 |
| Length | 64 | 64 | 64 | 64 |
| Height | 30 | 30 | 30 | 52 |
| Weight | 295 | 295 | 295 | 270 |