

Mobile Exam Tables

Transport, Lift, Exam

Veterinary Mobile Exam Table

This exam table conveniently doubles as an exam and transport table. This all purpose mobile lift table gives you and your staff the ability to lift, move and examine your patient. Our mobile platform enables this high quality exam table to function safely as an emergency transport.

FEATURES


- Exam top is fabricated with 304 grade, 16 gauge stainless steel
- Top measures 22"W x 58"L with raised edge to prevent dripping
- Top has four tie-down restraint cleats and corner bumper guards
- The base is double powder coated and will stand up to cleaning and disinfecting solutions
- Foot operated pedals on both sides, raises from 17" to 44" high
- Smooth lowering and lifting function capacity up to 300 lbs
- A strong rechargeable battery pack maintains a 24 volt actuator system
- Four corrosion-resistant locking castors
- Optional nonstick table topper



POWER LIFT WITH FOOT CONTROL



EASY ACCESS: IMPRESSIVE LOW HEIGHT ADJUSTMENT



OPTIONAL: NONSTICK TABLE TOPPER

Stainless Mobile Exam Table Specifications

Vetline Series Comparison Chart

VetLine

Mobile Exam & Lift Table

Base Info

Lowers to: 17"

Raises to: 44"

Capacity: 300 lbs.

Power: Battery

Material: 7 ga.

Finish: Powder Coated

Top Specs

Width: 22"

Length: 58"

Material: 16 ga ss

Type: Fixed

Shape: Rectangle

Product Dims & Weight

Width: 22"

Length: 58"

Weight: 203 lbs.

Shipping Dims & Weights

Width: 40"

Length: 64"

Height: 30"

Weight: 295 lbs.

PETLIFT LINE:	VetLine Series	VetLine Series	VetLine Series	VetLine Series
PRODUCT:	The Mobile Exam	Economy Exam HE	LowMax Exam ET	Classic Exam
BASE INFO				
Lowers to	17	19	13 1/2	27
Raises to	44	44	41	37
Weight Capacity	300	250	300	200
Power	Battery	Hydraulic	Electric	Hydraulic
BASE SPECS				
Base Material	7 ga.	7 ga.	7 ga.	7 ga.
Base Finish	Powder Coated	Powder Coated	Powder Coated	Laminate
Standard Color	Black or Silver Vein	Black or Silver Vein	Black or Silver Vein	Natural Almond
TOP SPECS				
Left Hand(LH)/Right Hand(RH)	N/A	N/A	N/A	N/A
Width	22	22	22	22
Length	58	44	58	44
Shape	Rectangle	Rectangle	Rectangle	Rectangle
Top Material	16 ga ss	16 ga ss	16 ga ss	16 ga ss
Top Type	Fixed	Fixed	Fixed	Fixed
OVERALL DIMENSIONS				
Width	22	22	22	22
Length	58	44	58	40
CABINET SPECS				
No. of Drawers	N/A	N/A	N/A	1
No. of Doors	N/A	N/A	N/A	N/A
No. of Adjust Shelves	N/A	N/A	N/A	N/A
Miscellaneous	N/A	N/A	N/A	Rotating Base
ACCESSORIES INCLUDED				
Duplex Electrical Outlet on Base	No	No	Yes	No
Corner Bump Guards A19-BG18	Yes	No	No	No
Tie Downs A18-TD53	Yes	No	No	No
Adjust Floor Levelers	No	Yes	Yes	No
Locking Casters A20-100-45	Yes	No	No	No
Digital Scale 400 lbs. A16-PS500	No	No	No	No
ACCESSORIES AVAILABLE				
Duplex Electrical Outlet on Base	No	No	Incl.	No
Corner Bump Guards A19-BG18	Incl.	No	No	Yes
Tie Downs A18-TD53	Incl.	Yes	Yes	Yes
Adjust Floor Levelers A22-FL54	No	Incl.	Incl.	No
Locking Casters A20-100-45	Incl.	Yes	No	No
Waste Push Thru Flap A-20-100-45	No	No	No	No
Table Mats	Yes	Yes	Yes	Yes
COLOR OPTIONS				
Silver Vein	TEZ-9010BT-SV	TEZ-9025H-SV	TEZ-9095E-SV	N/A
Black	TEZ-9010BT-BLK	TEZ-9025H-BLK	TEZ-9095E-BLK	TBP-1010H-BLK
Natural Almond	N/A	N/A	N/A	TBP-1010H-NA
Designer White	N/A	N/A	N/A	TBP-1010H-DW
Fashion Gray	N/A	N/A	N/A	TBP-1010H-FG
Wood Grain: Kensington Maple	N/A	N/A	N/A	TBP-1010H-WGKM
Wood Grain: Wild Cherry	N/A	N/A	N/A	TBP-1010H-WGWC
SHIPPING DIMS & WEIGHTS				
Width	40	40	40	40
Length	64	64	64	64
Height	30	30	30	52
Weight	295	295	295	270