

# TRANSPORT GURNEYS

**INFINIUM  VET**

## Mobile Gurney Supply Transport Table with Bumpers


This gurney/mobile table is constructed of heavy-duty 16 gauge 304 stainless steel and includes a fully welded under shelf. The legs are constructed of 1-1/2" square stainless steel tubing with rolling bumpers/ corner bumpers and heavy-duty 5" diameter swivel casters (two with locking brakes). 10" under shelf height.

- 48"L x 24"D x 36"H
- 60"L x 24"D x 36" H
- 48"L x 24"D x 36" H ( MRI room compatible )

## Mobile Gurney with Anesthesia Pole Insert

- 48"L x 24"D x 36"H
- 54"L x 26"D x 36"H ( with pole insert )

### Optional E-Tank Bracket



**INFINIUM<sup>®</sup>**

Infinium Veterinary  
12151 62nd St North #5  
Largo, FL 33773 • USA

 1(866) 846-7517  
 [infiniumvet.com](http://infiniumvet.com)  
 [sales@infiniumvet.com](mailto:sales@infiniumvet.com)