

# OVER INSTRUMENT TABLES

## INFINIUM VET PREMIUM SERIES

This table is used to keep surgical instruments, supplies, and monitoring equipment close and at hand during most procedures. The over-instrument table features a telescoping base for versatility. with a crank handle that can be used on either side of table. The top has a slight beveled marine edge to prevent spillage or roll off of supplies as well as a 1" splash on the rear of table rear to prevent roll off of supplies.

- 39" to 61" telescoping base
- 40" wide by 23" deep
- 1" splash and beveled marine edge
- 304 type 16 grade stainless steel
- (4) 3" nonmarking casters (2) with brakes

## Standard SERIES

All stainless steel flat top that lowers and raises with removable easy turn crank handle cranks from both sides of the table.

- 38" to 59" telescoping base
- 40" wide by 23" deep
- 304 type 16 grade stainless steel

