

D-Series Data Logging Blanket & Fluid Warming Cabinets




Model # DWC243074T

(T = Chambers in Thirds)


DIMENSIONS

SHOWN WITH STANDARD 4" BASE




DOOR REMOVED FOR CLARITY

FRONT VIEW



SIDE VIEW



SPECIFICATIONS

Overall Size

26.5"D x 30"W x 74.5"H
(673.1mm x 762mm x 1892mm)
17.21 Cubic Feet Storage
3 Adjustable Shelves

Upper Chamber I.D.

23"D x 26"W x 15.25"H
(584.2mm x 660.4mm x 387.3mm)
5.28 Cubic Feet Storage
Temperature Range 90°-160° F (Tolerance +/- 1° F)

Lower Chamber I.D.

23"D x 26"W x 34.5"H
(584.2mm x 660.4mm x 876.3mm)
11.93 Cubic Feet Storage
Temperature Range 90°-160° F (Tolerance +/- 3° F)

Electrical

120V, 12.5 Amp, 50/60 Hz, .90 kWh (Avg.),
3071 BTU/hr (Avg.)

Door Types

Stainless Steel or Glass

MAC Medical, Inc.
Publication No. SPEC-16 Rev D
January 2023

Information regarding this product is
subject to change without prior notice.