


D-Series Data Logging Blanket & Fluid Warming Cabinets



Model # TWC183078


DIMENSIONS

SHOWN WITH
STANDARD 4" BASE




DOOR REMOVED
FOR CLARITY

FRONT VIEW



SIDE VIEW



SPECIFICATIONS

Overall Size

20.5"D x 30"W x 78.75"H
(520.7mm x 762mm x 2000mm)
10.36 Cubic Feet Storage
1 Adjustable Shelf

Upper Chamber I.D.

17"D x 26"W x 11.25"H
(431.8mm x 660.4mm x 285.7mm)
2.88 Cubic Feet Storage
Temperature Range 90°-110° F (Tolerance +/- 1° F)

Intermediate Chamber I.D.

17"D x 26"W x 10.75"H
(431.8mm x 660.4mm x 273mm)
2.75 Cubic Feet Storage
Temperature Range 90°-150° F (Tolerance +/- 1° F)

Lower Chamber I.D.

17"D x 26"W x 18.5"H
(431.8mm x 660.4mm x 469.9mm)
4.73 Cubic Feet Storage
Temperature Range 90°-160° F (Tolerance +/- 1° F)

Electrical

120V, 12.9 Amp, 50/60 Hz, .90 kWh (Avg.),
3071 BTU/hr (Avg.)

Door Types

Stainless Steel or Glass

MAC Medical, Inc.
Publication No. SPEC-19 Rev B
January 2023

Information regarding this product is
subject to change without prior notice.