


D-Series Data Logging Blanket & Fluid Warming Cabinets



Model # TWC243078


DIMENSIONS

SHOWN WITH STANDARD 4" BASE




DOOR REMOVED FOR CLARITY

FRONT VIEW



SIDE VIEW



SPECIFICATIONS

Overall Size

26.5"D x 30"W x 78.75"H
(673.1mm x 762mm x 2000mm)
14.01 Cubic Feet Storage
1 Adjustable Shelf

Upper Chamber I.D.

23"D x 26"W x 11.25"H
(584.2mm x 660.4mm x 285.7mm)
3.89 Cubic Feet Storage
Temperature Range 90°-110° F (Tolerance +/- 1° F)

Intermediate Chamber I.D.

23"D x 26"W x 10.75"H
(584.2mm x 660.4mm x 273mm)
3.72 Cubic Feet Storage
Temperature Range 90°-150° F (Tolerance +/- 1° F)

Lower Chamber I.D.

23"D x 26"W x 18.5"H
(584.2mm x 660.4mm x 469.9mm)
6.4 Cubic Feet Storage
Temperature Range 90°-160° F (Tolerance +/- 1° F)

Electrical

120V, 12.9 Amp, 50/60 Hz, .90 kWh (Avg.),
3071 BTU/hr (Avg.)

Door Types

Stainless Steel or Glass

MAC Medical, Inc.
Publication No. SPEC-20 Rev D
January 2023

Information regarding this product is
subject to change without prior notice.