


FLUOROKENNEL

INFINIUM
VETERINARY HEALTH



**A Dramatically Safer and More
Efficient Method for Digestive Imaging**

FluoroKennel

Digestive Imaging Technology



A precision engineered kennel providing a natural alternative to recumbency and restraint techniques during **fluoroscopic aerodigestive and swallowing disorder studies** and imaging. Free-standing, self-feeding, less handling, and less anxiety sets FluoroKennel apart from other assessment and C-ARM assistance technology.

- ✔ Permits assessment of natural feeding behavior.
- ✔ Provides safe, low stress space sized matched for animal comfort.
- ✔ Reduces mess and minimizes movement artifact to permit quality videofluoroscopic & spot film imaging.
- ✔ Decreases the number of essential personnel by eliminating the need for physical restraint.
- ✔ Restraint-free testing prevents radiation exposure to personnel while minimizing aspiration risk associated with syringe and force feeding methods.

**SAFE, EFFICIENT, & EFFECTIVE
DIGESTIVE STUDIES**



INDUSTRIAL STRENGTH



RADIOTRANSLUCENT




Natural Self Feeding & Free-Standing


- ✓ Transparent behavior-based design.
- ✓ Trapezoidal shape with wide entry door that naturally guides the animal inside.
- ✓ Food bowl partition controls access to food / liquid items during testing.
- ✓ Bowl rings are adjustable to accommodate animals of different sizes and feeding position.


PERSONALIZED FIT

- ✓ Three standard designs to accommodate various sized companion animals.
- ✓ Removable partitions easily adjust the interior length of the kennel for added fit.




3 Partition Slots






Front & Side Entry
with Sliding Access Door



FluoroKennel

Digestive Imaging Technology



3 FluoroKennel™ sizes accommodating 9 animal body length dimensions

Dimensions	SMALL			MEDIUM			LARGE		
<i>Wither Height</i>	12"	12"	12"	20"	20"	20"	28"	28"	28"
<i>Body Length</i>	12"	16"	22"	20"	26"	32"	30"	36"	44"

All dimensions are in inches. Insert appropriate FluoroKennel™ partition to reduce internal length.

Unit Dimensions

	Length (inches)	Width (inches)		Height (inches)	Weight (pounds)	
		Front	Back		Kennel	with partitions
Small	24	6	12-13	16	31	40
Med	36	7	18-19	24	61	77
Large	48	8	23-25	32	100	126



INFINIUM VETERINARY HEALTH

Infinium Veterinary
12151 62nd St North #5
Largo, FL 33773 • USA



📞 1(866) 846-7517
🌐 infiniumvet.com
✉️ sales@infiniumvet.com