

STAINLESS STEEL WET TABLES

INFINIUM  VET

All Stainless Steel Cabinet Style Wet Table with Knee Space

Constructed of 16 ga. Type 304 stainless steel
Constructed of heavy-duty 16 gauge 304 stainless steel.
The wet table insert is 6" deep and includes S/S prep rack.
Features hinged cabinet door and three drawers with stainless steel pans. Also available with two cabinet doors in lieu of drawers. The drain is located at the center, by the faucet.

Model Options

- 000.933.0054** With one cabinet door and three drawers
Dimensions: 60"L x 24.5"D x 36"H
- 000.933.0055** With two cabinet doors (no drawers)
Dimensions: 60"L x 24.5"D x 36"H
- 000.933.0056** With one cabinet door and three drawers
Dimensions: 48"L x 24.5"D x 36"H
- 000.933.0057** With two cabinet doors (no drawers)
Dimensions: 48"L x 24.5"D x 36"H

Please select **RIGHT** or **LEFT** knee space configuration.

Also ssk about plumbing: Faucets, sprayers, inserts, etc.



INFINIUM[®]

Infinium Veterinary
12151 62nd St North #5
Largo, FL 33773 • USA

 1(866) 846-7517
 infiniumvet.com
 sales@infiniumvet.com