

Surgical & Exam

Panno-med PURE Treatment & Exam Table Platform

Engineered for innovation, durability, and flexibility.

Panno-med PURE - Veterinary Adjustable Universal Lift Table is ideal for a variety of surgery, exam, and treatment applications. Its refined stainless steel materials and simplified design make it versatile, strong, and exceptionally easy to clean. It comes with a variety of tops and can be free standing or wall mounted. It comes with a powerful electric lift and tilt mechanism and can be adjusted from nearly floor level to over 44".

- Versatile premium surgery and treatment table co-developed with well-known respected experts in veterinary medicine
- High strength powder coated European steel frame for maximum stability
- Very easy to clean, tightly sealed design with minimal parts, and protected lip gaskets
- Electrical height adjustment by footswitch
- Minimum height of 4.3 in resp. 7.1 in (incl. stand)
- Maximum height of 44.1 in resp. 46.9 in (incl. stand)
- Table top options: Fixed top or tiltable up to 12° with one-hand gas spring operation
- Top plate options: Stainless steel or synthetic top with edge and drain (51.2 x 23.6 inches) as well as 14 tie-down hooks
- Free-standing or highly secure wall-mount model options
- Modular build for future adjustments and versatility
- Lifting capacity up to 264.6 lbs with uniform load distribution



Panno-Med PURE Optional Configurations

Optional Weight Scale Features

- 0-220 lbs in 0.1 lbs steps
- Animal weighing function (averaging method)
- Battery-operated (220 hours working hours)
- Automatic shutdown after 3 minutes
- Delivery without batteries (6 pieces AA)

Free Standing Version Top Options

- Fixed stainless steel top
- Tilting stainless steel top
- Fixed synthetic top
- Tilting synthetic top
- Fixed stainless steel top with integrated scale
- Fixed synthetic top with integrated scale



Wall-Mounted Version Top Options

- Fixed stainless steel top
- Tilting stainless steel top
- Fixed synthetic top
- Tilting synthetic top
- Fixed stainless steel top with integrated scale
- Fixed synthetic top with integrated scale

